

H2AFJ-H2A histone family, member J Gene

Reference: PTXBC003602

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2AJ; histone H2A.J; H2a/j; buforin I; H2A histone family member J		

For research use only.

GENE DESCRIPTION	H2A histone family, member J
GENE ID	55766
ORIGIN SPECIES	Human
ACCESSIONS	BC003602
VECTOR	pDONR223
NCBI SYMBOL	H2AFJ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccggtcgggaaacaggcggcgaagtgcgagcaaggccaatccgctcctcccgcgccgctgcagttccggtggccgagtgacagactgctgcgcaagggaactacgggagcgagtgggcgccggcgccggtgtacctggcgcggtgttgagta
ccttacggcgagatcctggagctggctggcaacgcccgcgtgacaacaagaagaccaggataatccccgccacctgcagctcgccatccgcaacgacgagaggttaacaagctgctgggcaagtgcacatgctcaggcgccgctcctgccaacatccaggccgtgctgct
gccaagaagacggagagtcagaagacgaagagcaaatga