



Reference: PTXBC003635

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MPSL1; PUMP-1; matrilysin; matrin; matrix metalloproteinase 7 (matrilysin, uterine); matrix metalloproteinase 7 (matrilysin, uterine); matrix metalloproteinase-7; pump-1 protease; uterine matrilysin; uterine metalloproteinase; matrix metalloproteinase 7		

For research use only.

GENE DESCRIPTION	matrix metalloproteinase 7 (matrilysin, uterine)
GENE ID	4316
ORIGIN SPECIES	Human
ACCESSIONS	BC003635
VECTOR	pDONR223
NCBI SYMBOL	MMP7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]