

CDK4-cyclin-dependent kinase 4 Gene

Reference: PTXBC003644

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMM3; PSK-J3; cyclin-dependent kinase 4; cell division protein kinase 4; cyclin dependent kinase 4		

For research use only.

GENE DESCRIPTION	cyclin-dependent kinase 4
GENE ID	1019
ORIGIN SPECIES	Human
ACCESSIONS	BC003644
VECTOR	pDONR223
NCBI SYMBOL	CDK4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctacctctgatatgagccagtggtgaaattggtgctggtgacatgggacagtgtaacagccggtgacccacagtgccactttgtgcccctcaagagtgtgagagtcaccaatggaggaggaggtggaggaggcctccatcagcacagttcgtgaggtggttactga
ggcgactggaggcttttgagcatcccaatgttgccggtgatggacgtctgtgccacatcccgaactgaccgggagatcaaggtaacctgtgtttgagcatgtagaccaggacctaaggacatctggacaaggcaccaccaggttgcagccgaacgatcaaggtctga
tgcgccagtttctaagaggcctagatttctcatgccaattgcatcgttcaccgagatctgaagccagagaacattctgtgacaagtgtgtgaacagtaagctggtgactttggcctggccagaatctacagctaccagatggcacttacaccgtgtgttacactctggtaccgagct
cccgaagtcttctgcagtcacatgatcaacacctgtggacatgtggagtgttggctgtatctttgcagagatgtttcgtgaaagcctctctctgtggaactctgaagccgaccagttggcgaatctttgacctgattgggctgctcagaggatgactggcctcagatgtatccctgc
cccgtggagcctttccccagaggccccgccagtgagtcggtgtacctgagatggaggagtcgggagcacagctgctgtgaaatgctgacttttaaccacacaagcgaatctctgctttcgagctctgcagcactcttatcataaggtgaaggtatccggagtga