

PDXK-pyridoxal (pyridoxine, vitamin B6) kinase Gene

Reference: PTXBC003651

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf124; C21orf97; HEL-S-1a; PKH; PNK; PRED79; epididymis secretory sperm binding protein Li 1a; pyridoxamine kinase; pyridoxine kinase; vitamin B6 kinase; pyridoxal (pyridoxine, vitamin B6) kinase		

For research use only.

GENE DESCRIPTION	pyridoxal (pyridoxine, vitamin B6) kinase
GENE ID	8566
ORIGIN SPECIES	Human
ACCESSIONS	BC003651
VECTOR	pDONR223
NCBI SYMBOL	PDXK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccctccatctctgaccttgccctgggacaagctttgatggggggcccccagcttcaagctgtgtgtgggaacagcaccceaaatgccagctctctcttccatccaccagtatactgctggggccattctgtgtctttgtccaacagaaacccatttctgtgggatgcctcc
agtgccacagggccactcaccatgcatctctgtctgcccgtcagtgctgggacggagacagcaaggcgaagccagtgctgtggcagataggtgggtgggaacagagaggggagaatgccgtcctaagcttctgtctgggatccccacacacacctgggtactgctgggaacctg
tcctaagtaaaactatggacctgcctcgccaccggcctgcgaagccagcatctcctgtaagggtgatggaagcgcctttgtctcactttgagctgcaagctgggtcagcggctctgaagccctcagtgactttctaaccaagaccagccctggcaggaggaggggtgggtgcag
ggctggtgggacaaaagaggcctcagcaggcctggaagacccttcagtcacatccacagcgtgtcagcagctgggagaacctgtgtcaagctcgagccgtcataggtcccatgaggtgtctgaagcccttctgtgtga