

ATP5J2-ATP synthase, H⁺ transporting, mitochondrial F0 complex, subunit F2 Gene

Reference: PTXBC003678

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP5JL; ATP synthase subunit f, mitochondrial; ATP synthase f chain, mitochondrial; ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit f; F1F0-type ATPase subunit f; F1Fo-ATP synthase complex Fo membrane domain f subunit; F1Fo-ATPase synthas		

For research use only.

GENE DESCRIPTION	ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit F2
GENE ID	9551
ORIGIN SPECIES	Human
ACCESSIONS	BC003678
VECTOR	pDONR223
NCBI SYMBOL	ATP5J2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtcagttggtagtgcggcccgaccagtgaggacaagaactctggaggtcaactctggagctgccaagctggatcttgatcgaggacttcagtcctagtgccatttcggagcgtttcaagaggttactaccggtactacaacaagtacatcaatgtgaagaagggg
agcatctcggggattaccatgggtgctgcatgctacgtgctcttttagctactcttttctacaagcatctcaagcagcagcggtctccgcaataccactga