

MED6-mediator complex subunit 6 Gene

Reference: PTXBC004106

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARC33; NY-REN-28; mediator of RNA polymerase II transcription subunit 6; CTD-2540L5.5; activator-recruited cofactor 33 kDa component; renal carcinoma antigen NY-REN-28; mediator complex subunit 6		

For research use only.

GENE DESCRIPTION	mediator complex subunit 6
GENE ID	10001
ORIGIN SPECIES	Human
ACCESSIONS	BC004106
VECTOR	pDONR223
NCBI SYMBOL	MED6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggtggatccgagacaatctgctgggaattcttgggttgacagctcttgatccctatttgaacagtggttagtgcctggattactttcagaaagaagtaaccttttatgacagaacatgtaataatgaagtgggtcaaaatgcagaggctaaccattagaacacttgaatcagat
gggtggaaatcagatcacatctttgcatgctcaagagcccattctttcatcattcggagcaacagcggcagctccctgcccaagttatccactagctgattactatattgctggagtgatctatcaggcaccagactgggatcagttataaactctagagtgcttactgagtcagtgat
tcagtcagcttttgatgaagctatgcatcactgctgataatccttccaaaggggtattggtggcacttcaagatcatgaagagcaagataaagtcagacctaaagccaaaaggaaagaacaaagctctattttcagagacaacgtgtggatgcttcttttagacctcagacaaaattt
ccacccaaatttgtagcagctaaagcctggagaaaaagcctgttccagtgatcaacaagaagaggcagaacctataccagaactgtaaaacctgaggagaaggagaccacaagaatgtacaacagacagtgagtgctaaaggccccctgaaaaacggatgagacttcagtga