

C22orf13-chromosome 22 open reading frame 13 Gene

Reference: PTXBC004144

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C22orf13; LLN4; protein GUCD1; CG13760 gene product [Drosophila melanogaster] homolog; guanylyl cyclase domain-containing protein 1; guanylyl cyclase domain containing 1		

For research use only.

GENE DESCRIPTION	chromosome 22 open reading frame 13
GENE ID	83606
ORIGIN SPECIES	Human
ACCESSIONS	BC004144
VECTOR	pDONR223
NCBI SYMBOL	C22orf13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaccatgccggcctaattttttgtatttttagtagaccaggttcgccatgttgccaggctgtctcgaacccctgagctcaggtgatccacctgctcggcctccaaagtgggaatgggtcattttcaaaggggcaacagtataccatccctgcactcctgggtcctgtcttc
caatgcccaaggcatggcatgggtgagctgtgccccgccattgtttggagcacagatga