

LYZ-lysozyme (renal amyloidosis) Gene

Reference: PTXBC004147

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LYZF1; LZM; lysozyme C; 1,4-beta-N-acetylmuramidase C; c-type lysozyme; lysozyme F1		

For research use only.

GENE DESCRIPTION	lysozyme (renal amyloidosis)
GENE ID	4069
ORIGIN SPECIES	Human
ACCESSIONS	BC004147
VECTOR	pDONR223
NCBI SYMBOL	LYZ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggctctcattgttctgggcttgcctcttctgttacgggtcagggcaaggctttgaaaggtgtgagttggccagaactctgaaaagattgggaatggatggctacaggggaatcagcctagcaaaactggatgtgttggccaaatgggagagtggttacacacagagctacaa
actacaatgctggagacagaagcactgattatgggatattcagatcaatagccgctactggtgtaatgatggcaaaacccagagcagtaaatgcctgtcatcttcctgcagtgctttgctgcaagataacatcgtgatgctgtagcttgtgcaagagggtgtccgtgatccacaaggc
attagagcatgggtggcatggagaaatcgtgtcaaaacagagatgtccgtcagtatgtcaaggttgggagtgtaa