

B9D2-B9 protein domain 2 Gene

Reference: PTXBC004157

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ICIS-1; MKS10; MKSR2; B9 domain-containing protein 2; MKS1-related protein 2; involved in cilia stability-1; B9 protein domain 2		

For research use only.

GENE DESCRIPTION	B9 protein domain 2
GENE ID	80776
ORIGIN SPECIES	Human
ACCESSIONS	BC004157
VECTOR	pDONR223
NCBI SYMBOL	B9D2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgaggtgcacgtgatcgggcagatcatggggccagcgggtttctcggaaagtagcctctctgcaagtggggcattcacacagggcgcatggaagctcctgtcagcgctgcgggagggccaaacgcaagtggacaccccgagataggggacatggcttactgtgccac
cccatcgacctgcacttcgcccaaaaggtctcaaggctggccccggctccattccagggtgtgtgccaggacagctttggccgctgccagcttgcaggctatggaattttgcatgtgccagtagccggggcaccaccagctggcctgccccacgtggcggccccctgggcagttgg
cgagaacagttggcacgggctttcgtgggtgtgggcccagctgctgcatggggacaccatctacagtgggggcaccgctatcgctgcacacagctgctgggtggcacctggagatcggcctgctcgcgaacttcgaccgctacggcgtggagtgctga