

DCUN1D5-DCN1, defective in cullin neddylation 1, domain containing 5 (S. cerevisiae)
Gene

Reference: PTXBC004169

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SCCRO5; DCN1-like protein 5; DCN1, defective in cullin neddylation 1, domain containing 5; DCUN1 domain-containing protein 5; defective in cullin neddylation protein 1-like protein 5; defective in cullin neddylation 1 domain containing 5		

For research use only.

GENE DESCRIPTION	DCN1, defective in cullin neddylation 1, domain containing 5 (S. cerevisiae)
GENE ID	84259
ORIGIN SPECIES	Human
ACCESSIONS	BC004169
VECTOR	pDONR223
NCBI SYMBOL	DCUN1D5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggtgaagaagaagaaaaatccctgggtggcagcagcagtagcggagacggaggcccaaaaagtgtaaaatccagctattgcagatcccaacccctgctagactaataagtggagaggaaacattttcaagcaagaagtcctggctgtttatgaatgcaggtc
ctgatgaagttgtagggcagaaggaatgaaaaattgtgaagacattgtgttgacacgtgaaatattattattgttttagcgtggaattggagcctgaaagcatgggatttttaccaaggagaatggttaagggaatgacttcattacagtgtgactgcacagaaaagtacaaaa
caaatgtacttttgcgtcacagttgaatgataattcgtcatttaagaatactacagatagccttgattttgcaaggataaagatcagagaagccttgatattgatactgctaaatctatgtagctctctctgtgggagacatggccactgtttcagttattaccagtacctggagcaatcaa
agtatcgtgttatgaacaaatcaatggtacaatgtattagaattcagcagaacagtcacgtcgtgatcttagtaactatgatgaagatggtgcttgccctgttcttcttgatgaattgttgagtgccaaaaagtcctgcagacatcatag