



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SSH3-slingshot homolog 3 (Drosophila) Gene

Reference: PTXBC004176

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SSH3L; protein phosphatase Slingshot homolog 3; SSH-3L; SSH-like protein 3; hSSH-3L; slingshot 3; slingshot homolog 3; slingshot protein phosphatase 3		

For research use only.

GENE DESCRIPTION	slingshot homolog 3 (Drosophila)
GENE ID	54961
ORIGIN SPECIES	Human
ACCESSIONS	BC004176
VECTOR	pDONR223
NCBI SYMBOL	SSH3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr