



Antibodies and Proteins

Espace Européen de l'Entreprise

67300 Schiltigheim, France

contact@proteogenix.fr

www.ProteoGenix.science

Reference: PTXBC004204

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNA-binding protein 42; RNA binding motif protein 42		
For research use only.			

GENE DESCRIPTION	RNA binding motif protein 42
GENE ID	79171
ORIGIN SPECIES	Human
ACCESSIONS	BC004204
VECTOR	pDONR223
NCBI SYMBOL	RBM42
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]