

SIRT6-sirtuin (silent mating type information regulation 2 homolog) 6 (S. cerevisiae)
Gene

Reference: PTXBC004218

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SIR2L6; NAD-dependent protein deacetylase sirtuin-6; SIR2-like protein 6; regulatory protein SIR2 homolog 6; sir2-related protein type 6; sirtuin type 6; sirtuin 6		

For research use only.

GENE DESCRIPTION	sirtuin (silent mating type information regulation 2 homolog) 6 (S. cerevisiae)
GENE ID	51548
ORIGIN SPECIES	Human
ACCESSIONS	BC004218
VECTOR	pDONR223
NCBI SYMBOL	SIRT6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagcagggctctggcccccaagttcgacaccactttgagagcgcgccgcccacgacacccacatggcgctgggtgcagctggagcggctgggctctccgcttctggtcagccagaacgtggaggggtccatgtgcgtcaggctcccaggagacaaactggcaga
gctccacgggaacatgtttgtggaagaatgtgccaagtgaagacgcagtagctccgagacacagtcgtgggcaccatgggctgaaggccacgggcccgtctgcacgtgggtaaggcaaggggctgcgagcctcaggacgcccagctgtccatcacgtgggtacatcgt
gcagatccggccagcgggaacctgcccgtggtacccaagcggggagggcgctggtcatcgtcaacctgcagcccaccaagcacgaccgtatgctgacctccgcatccatggctacgttgacgaggtcatgacccggctcatgaagacctggggctggagatccccgct
gggacggccccgtgtgctgagagggcgctgccacctgccccgcccccccccaagctggagcccaaggagggaatctccacccggatcaacggctctatccccggcgcccccaaggaggagccctgcgcccagcacaacggctcagagcccgccagccccaaacg
ggagcggccaccagccctgccccccacagacccccaaaagggtgaaggccaaggcggtccccagctga