

LAMA4-laminin, alpha 4 Gene

Reference: PTXBC004241

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LAMA4*-1; CMD1JJ; LAMA3; laminin subunit alpha-4; laminin alpha 4 chain; laminin, alpha 4; laminin subunit alpha 4		

For research use only.

GENE DESCRIPTION	laminin, alpha 4
GENE ID	3910
ORIGIN SPECIES	Human
ACCESSIONS	BC004241
VECTOR	pDONR223
NCBI SYMBOL	LAMA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctttgagctcagcctggcgtcgttctgctctgtggtcctctgagcgcgtgctgtcccgccgcgtccgggacgacaacgttttctttgacattgaaggagctcagcgttggcaggcaagacccgctgagacgagcgaaccccgctggctctgggacgcctg
ccgcctgcggccgaggtacagtgccctgccattgccaccctgctggggcactgcgccccgcgggctgtgcccacactcgtccttctctctcctcgcctcttctctccccagtgcttgagagtttcacctgggctaggtcagttcggaaacttgaataaagagtttctttgttaa