



Reference: PTXBC004248

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OP-1; bone morphogenetic protein 7; osteogenic protein 1		
For research use only.			

GENE DESCRIPTION	bone morphogenetic protein 7
GENE ID	655
ORIGIN SPECIES	Human
ACCESSIONS	BC004248
VECTOR	pDONR223
NCBI SYMBOL	BMP7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]