



Reference: PTXBC004257

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ret proto-oncogene; ret proto-oncogene (multiple endocrine neoplasia and medullary thyroid carcinoma 1, Hirschsprung disease); proto-oncogene c-Ret; RET transforming sequence; RET receptor tyrosine kinase; CUX1/RET fusion; proto-oncogene tyrosine-protein		

For research use only.

GENE DESCRIPTION	ret proto-oncogene
GENE ID	5979
ORIGIN SPECIES	Human
ACCESSIONS	BC004257
VECTOR	pDONR223
NCBI SYMBOL	RET
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]