

**PHF15-PHD finger protein 15 Gene**

Reference: PTXBC004292

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PHF15; JADE-2; protein Jade-2; PHD finger protein 15; jade family PHD finger protein 2; jade family PHD finger 2		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	PHD finger protein 15
<b>GENE ID</b>	23338
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC004292
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	PHF15
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atggccatgagcagtgagccactgaacaatggcaccgcgaggaccctgctccaggctgctgctcagaggaactgctcaggacgaggagacactgctcagcttcatcgggaccctcgtcgacctggtgaccctgctaggaagcccaggccgaccgcctgctgcca  
agaagaacaccaccaccaccaccgagacgggctggttcacggagcactccagacaagcccccaagaagacctggggccaggatcagcagtggaaggggtcaagggccacctacaggaagccaccacgtcggacatctctcacttccgtccagccctgag  
ccggggactgtccatctagccaccctgaaagccccccgactggccctgagaccgggacgaggcagcctcagtagctgctgactcagatgtccaagtgcctggccctgagcaagccctaaagccttgggcccgtccggccaccggagcaaggtaccggaga  
ttgccgggtgccagcctgatctgggatgggaccacctcagctgtggctgagaggccaagatcagcctgcatfttgactgagactgatgctacttctctgatgggagatgagcactcagatgtagagccgaggacggtgggtgcagcgggtccccgggagcagg  
gcagaggaggtgtccgatgggctactggctcctaa