

RPL3-ribosomal protein L3 Gene

Reference: PTXBC004323

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ASC-1; TARBP-B; 60S ribosomal protein L3; HIV-1 TAR RNA-binding protein B; ribosomal protein L3		

For research use only.

GENE DESCRIPTION	ribosomal protein L3
GENE ID	6122
ORIGIN SPECIES	Human
ACCESSIONS	BC004323
VECTOR	pDONR223
NCBI SYMBOL	RPL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgcctgtctcctctgcgcagaagaaggccacctgatggagatcaggtgaacggaggcactgtgcccagaagctggactggcccgagagggcttgagcagcaggtacgtgtgaaccaagtgtttgggcaggatgatgatcgcgtcatcggggtgaccaaggcga
aggctacaagggtcaccagtcgttggcacaccaagaagctgccccgaagaccaccgaggcctgcgaaggtggcctgtatfggggcatggcatcctctctgtgtagccttctctgtggcacgcgtgggcagaaggctaccatcacgcgactgagatcaacaagaattata
agattggccagggtacattatcaaggacggcaagctgatcaagaacaatgcctccactgactatgacatctgacaagagcatcaaccctctgggtggctttgtccactatggtgaagtgaacaaatgactttgtcatctgaaaggctgtgtgtgggaaccaagaagcgggtgctcacc
ctccgcaagtcctgtgtgagacgaagcggcggtctgtggaagattgaccttaagttcattgacaccacctcaagtttggccatggccgttcagaccatggaggagaagaagcattcatgggaccactgaagaagaccgaattgcaaggaagaaggagcttaa