

C12orf39-chromosome 12 open reading frame 39 Gene

Reference: PTXBC004336

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf39; spexin; NPQ; neuropeptide Q; spexin hormone		
For research use only.			

GENE DESCRIPTION	chromosome 12 open reading frame 39
GENE ID	80763
ORIGIN SPECIES	Human
ACCESSIONS	BC004336
VECTOR	pDONR223
NCBI SYMBOL	C12orf39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagggaactcagaagctctggcagcaacaacctggctcttttctctgtgtttttctctgggaaactccagctgcgtccgcagagactgttgagagaaggaaactggactctcaagctatgctctacctgaaaggggcacagggtcgccgttcatctccgaccagagccggagaaa
ggacctctccgaccggccactgccggaagacgaagcccaatcccaactactaactattccggaggcagcaaccatcttactgctgctctcagaaatcaccagaagatgaagaaaaaactttgatcaaacagattcttggaagacagtctgcttaactggtga