

C21orf122-chromosome 21 open reading frame 122 Gene

Reference: PTXBC004343

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MGC10960; PRED57; PRED90; chromosome 21 open reading frame 122		
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 122
GENE ID	728039
ORIGIN SPECIES	Human
ACCESSIONS	BC004343
VECTOR	pDONR223
NCBI SYMBOL	C21orf122
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcaggctgccgcttcagtcactgcctgggaaacaagcctgtctttggggccagccagctaacaagtgaagccctcaccatctgaaggcgaagcaccacaaaggggcaaccgtggacggctaa