

POLR3GL-polymerase (RNA) III (DNA directed) polypeptide G (32kD)-like Gene

Reference: PTXBC004355

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RPC32HOM; flj32422; DNA-directed RNA polymerase III subunit RPC7-like; DNA-directed RNA polymerase III subunit G-like; RNA polymerase III subunit C7-like; RPC32-like protein; alternative RNA polymerase III subunit 32; polymerase (RNA) III (DNA directed) p		

For research use only.

GENE DESCRIPTION	polymerase (RNA) III (DNA directed) polypeptide G (32kD)-like
GENE ID	84265
ORIGIN SPECIES	Human
ACCESSIONS	BC004355
VECTOR	pDONR223
NCBI SYMBOL	POLR3GL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagccgggtggggccgggtgctggccggggccagttgacctcaacgtggagggccgtgggcatgggaaaggggatgcttggcccccacccctgcagccttctccactcttccctccctggagttccgccagtagcttggccctcaggcgaggagggaatatgt
cctggcactgaagcaagagctacgaggagccatgaggcagctccctacttcatccggccagctgtcccaagagagatgtggagcgttattcagacaaatcatgatgtcaggtccgattgacaatgccatcgattggaaccctgattggcggtctacccgggagctaaagatcgg
agtgcggaagctacagaagggaacggattacaattctgctcccaagaggcccccctaagaccacagaagataaggagggaacatacagaactagagacccctggagaagaaggagaagaagtaactcagaggaggatgaggagaagaagaagaagaaggaaggaagg
aggagaagaagagtagatgaagaagaacatgaagggaactgattacatcatgcatatttgacaatggagaggacttgggtgacagtgatgacaatattggacgaggtatatactga