

TMEM176B-transmembrane protein 176B Gene

Reference: PTXBC004358

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LR8; MS4B2; transmembrane protein 176B; LR8-like protein		
For research use only.			

GENE DESCRIPTION	transmembrane protein 176B
GENE ID	28959
ORIGIN SPECIES	Human
ACCESSIONS	BC004358
VECTOR	pDONR223
NCBI SYMBOL	TMEM176B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgcaaaacacgggtgattggaatggagttgctatggcctctaggccatccagcccccacgtcaacgtccacatccaccaggagtcagctttgacacaactctgaaagctggaggttctctgaagaagttcttttaccctggggacactgtgcctccacagccaggattggt
atgagcagctggctctaggggtgactcagatattgctgggggttgtagttgtgtctggagtggtctcagcttggggccctggagctgtgctgctcagcagctgtgctccttgggcccgggtctgtggtgatcgagcagagctggggccattgtccatgagaagcaccgggcaaa
ctgtgctgctatataccagcctgctcaccctgacaggcttggctacagctatggctgtgtgtcctctgctggaatagcttcatctggcacaactgaacccttttatacatcgacactgtgtgtgatcgtcagacccctgtctcctaccactgggtacagatggatgcggcgaagtcaagagaa
ccaatggcagaaggaggagtgtagagcttacctgacagatgctgaggaagttgttcacagcaatccgtgcccctgttctggtgctgtgtgtctgaaggctattgtgtcttgggttccctgggagtaggtcttgcgaactgtgtgtgcccagagctccagccctgaatgaggaggaatcagag
aagaggctactgggggagaattcagtgcccccttgcctctaggagcagacctccactgcatgtctctgga