

TMEM39B-transmembrane protein 39B Gene

Reference: PTXBC004364

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 39B		
For research use only.			

GENE DESCRIPTION	transmembrane protein 39B
GENE ID	55116
ORIGIN SPECIES	Human
ACCESSIONS	BC004364
VECTOR	pDONR223
NCBI SYMBOL	TMEM39B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccacccatgctgctcctgtcaccagcctcatccgagtgagggtggagttcctcaagatggactcaactggcgcatgaaggaagtgctcgtcagctccatgctgagcgcctactatgtggccttggcctgtctgttgctgaagaacacacattactatgacaagcgtggctctgtgaactctctctgctggtgtccatcagcaccctccgtgatcctcatgcagcaccctgctgctgcccagctactgtgacctgctgcacaaggccgcccacatctgggctgttggcagaagtgaggaccagcgtgtgtctccaacgtgctgcagcaccctggactgaagaatgcatgtggccgcagggcgtgctggtgaagcacagcaagaagctctacaagccgttagggccactacaacgtggatccctctgacgtctccacttccgcttccattcttttcagcaaacccctgcggatcctcaacatcctcctgctgctggaggcgtgtgcatgtgtatcagctgtactccctaagtcctctgaaaagtgccaccagaccatctcgtggcctcatctctcagcaactactatgcctcttcaagctgtccgggaccgctgtgtattgggaaggcctactcactctgctagccccagagagacgtggaccaccttctctga