

MRPL49-mitochondrial ribosomal protein L49 Gene

Reference: PTXBC004378

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf4; L49mt; MRP-L49; NOF1; 39S ribosomal protein L49, mitochondrial; neighbor of FAU; next to FAU; mitochondrial ribosomal protein L49		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L49
GENE ID	740
ORIGIN SPECIES	Human
ACCESSIONS	BC004378
VECTOR	pDONR223
NCBI SYMBOL	MRPL49
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagctaccatgttccggctacgtcggggatggagaaccggtgtccagcggggctgcgggctacggctgttgagccagaccagggccctccagattaccaggtttgtgagctctgtgatgaatatcagtttggagcgcctgttacggctaccagatccagatccc
ccaaagcatgaacattatctacccctagtggctggcagctccagagacccccacccaacctgccttctttgtacgacgtctcggatgcacaacatccccgtctacaaggacatcacgcatggcaaccggcagatgactgtgatccggaaagtggagggaacatctggccctg
cagaaagacgtggaagattttctgagcccgctgctggggaagacacctgtcaccaggtcaatgaggtgacaggtaccctacggatcaagggtactttgaccaggagcttaaagcctggctcttggaagaagccttctga