

MGC10981-hypothetical protein MGC10981 Gene

Reference: PTXBC004397

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AFAP1-AS; AFAP1AS; AFAP1 antisense RNA 1		

For research use only.

GENE DESCRIPTION	hypothetical protein MGC10981
GENE ID	84740
ORIGIN SPECIES	Human
ACCESSIONS	BC004397
VECTOR	pDONR223
NCBI SYMBOL	MGC10981
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaaggacactgccagcctcttcattcccctgtgtgatgctccaactccaaagtggaaaccagggcaaatgaatgtcggcagctgtcggatggcagcctcgtcttcgacagctggtgtgggggtaccacctggggctagctggcctcatcttcgttag