

C9orf64-chromosome 9 open reading frame 64 Gene

Reference: PTXBC004407

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0553 protein C9orf64; chromosome 9 open reading frame 64		
For research use only.			

GENE DESCRIPTION	chromosome 9 open reading frame 64
GENE ID	84267
ORIGIN SPECIES	Human
ACCESSIONS	BC004407
VECTOR	pDONR223
NCBI SYMBOL	C9orf64
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacctgggtgaaagtcttcttctacagagatgtgactctgtttgaggggaaaaagagttcttttacaacgagcccaatcctgtagcagatacgtggagtgtattggaaggaaaaggagatggctgctcaaggacatccagatcaccatgtttgctgattatagattacctcag
gttctgtctcatctggagccctgaaatactctgatgacctactgaagaagcttctcaaggagaaatgctctcatatggagacaggcaagaggtggaaatcagaggggtgctcgtttggtgtgttgagctgatccgggattgtctctggagcttattgaacaaaagggtgaaaaacctaagg
agagatcaattccattctctggtattattactatgggactatgcccatgacataggggaagatatgaagggaattccgtttcatcgcatacgttgcatatattattga