

C16orf57-chromosome 16 open reading frame 57 Gene

Reference: PTXBC004415

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0406 protein C16orf57; C16orf57; HVSL1; Mpn1; hUsb1; U6 snRNA phosphodiesterase; HVSL motif containing 1; U six biogenesis 1; U6 snRNA biogenesis 1; mutated in poikiloderma with neutropenia protein 1; U6 snRNA biogenesis phosphodiesterase 1		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 57
GENE ID	79650
ORIGIN SPECIES	Human
ACCESSIONS	BC004415
VECTOR	pDONR223
NCBI SYMBOL	C16orf57
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcgcggcgccctggtgggtacagcagcagcggctccgaggtgatggtccgagcgggatgcggaccagccggggatgggagccaccgtcgtggccagagcccccctccaggcagagattccagctacagtggtgctgaacatgttccgggaccgagggag
gggctgaagatgacagcacaacacgggggacgggtgcaccttccccacgagcagggcaactgggcccacacgtctatgtaccatataagccaaggaggagttcctggtctctgtatgtgtgtgctgccccatgccagacatatgtccccggctgtaaggatgaagg
tgttccacctcagcctgtccagagtggtgtctgcgccaccactggatcctccctctgtgcaggctctgaaagcccgatgacctctccacagattcttcttactgcccaaccaggtaagattacaccaatcaagagaaaaccaggaccttattgggcttgaggtcacttcaggcctatg
ccagttcctggacctggttcagagtggaacagagtcagtgagggaattcaacctcaccatttclaccagatccttctccacctcagcctggcctggtgtgtgggtgatgcacgtctccagctggaggggagtcctgcaggaactacaggcaatcgtggatgggttgaagatgctg
agggtcgtcgtcgcgtgcacactgagcaagtcgctgcaagtcgtgggaacaagttctctcgtatgcctttgaagtga