



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FUZ-fuzzy homolog (Drosophila) Gene

Reference: PTXBC004445

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NTD; protein fuzzy homolog; fuzzy homolog; fuzzy planar cell polarity protein		
For research use only.			

GENE DESCRIPTION	fuzzy homolog (Drosophila)
GENE ID	80199
ORIGIN SPECIES	Human
ACCESSIONS	BC004445
VECTOR	pDONR223
NCBI SYMBOL	FUZ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]