

PYCARD-PYD and CARD domain containing Gene

Reference: PTXBC004470

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CARD5; TMS; TMS-1; TMS1; apoptosis-associated speck-like protein containing a CARD; caspase recruitment domain-containing protein 5; target of methylation-induced silencing 1; PYD and CARD domain containing		

For research use only.

GENE DESCRIPTION	PYD and CARD domain containing
GENE ID	29108
ORIGIN SPECIES	Human
ACCESSIONS	BC004470
VECTOR	pDONR223
NCBI SYMBOL	PYCARD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcgcgcgcgacgccatcctggatgcgctggagaacctgaccgccgaggagctcaagaagtcaagctgcaggcgccacgcaccagggtctggagccgcgccagctgggatccaggccccctcctcagtcggcagccaagccaggcctgcactttatagaccagcac
cgggctgcgttatcgcgagggtcacaacgttgagtggtgctggatgctctgtacgggaaggtcctgacgagatgagcagtaaccaggcagtgccggccgagccccaccaaccaagcaagatgcggaagctctcagtttcacaccagcctggaactggacctgcaaggacttgctct
ccaggccctaaggagtgccagtcctacctggtggaggacctggagcggagctga