

HYOU1-hypoxia up-regulated 1 Gene

Reference: PTXBC004560

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GRP-170; Grp170; HSP12A; ORP-150; ORP150; hypoxia up-regulated protein 1; 150 kDa oxygen-regulated protein; 170 kDa glucose-regulated protein; oxygen regulated protein (150kD); hypoxia up-regulated 1		

For research use only.

GENE DESCRIPTION	hypoxia up-regulated 1
GENE ID	10525
ORIGIN SPECIES	Human
ACCESSIONS	BC004560
VECTOR	pDONR223
NCBI SYMBOL	HYOU1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagacaaagttaggagcagagggccgaggagcagctgtgtgggccttgggtgctgtctctggcagacctgtggcactgagtgatacactggcagtgatgtctgtggacctgggcagtgagtcctgaagggtggccattgtcaaacctggagtgcccatggaaattgtctga
ataaggaaatctcgaggaaaacaccgtgatctgaccctgaagaaaaatgaagattcttggagacagtgacgaagcatggcgattagaatccaaaggctacgtacttaccagcacctctctggggaagcaggcagataacccccatgtagctctttaccaggcccgcttccc
ggagcacgagctgacttgcaccacagaggcagactgtgcacttcagatcagctcgagctgcagttctcaccctgtga