

MRPS21-mitochondrial ribosomal protein S21 Gene

Reference: PTXBC004566

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MDS016; MRP-S21; RPMS21; 28S ribosomal protein S21, mitochondrial; S21mt; mitochondrial 28S ribosomal protein S21; mitochondrial ribosomal protein S21		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S21
GENE ID	54460
ORIGIN SPECIES	Human
ACCESSIONS	BC004566
VECTOR	pDONR223
NCBI SYMBOL	MRPS21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaaacatctgaagttcatcgccaggactgtgatgtacaggaagggaacgtggaaagcgcatcacaggaccctaaacagaatcctcactatggatgggctcattgaggacattaagcatcggcggtattatgagaagccatgccgccggcgacagagggaagctatgaaaggt
gccggcggatctacaacatggaaatggctcgcaagatcaactcttgatgcgaaagaatcgggcagatccgtggcagggctgctga