



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
 contact@proteogenix.fr
 www.ProteoGenix.science

FREQ-frequenin homolog (Drosophila) Gene

Reference: PTXBC004856

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FREQ; FLUP; neuronal calcium sensor 1; frequenin homolog; frequenin-like protein; frequenin-like ubiquitous protein		

For research use only.

GENE DESCRIPTION	frequenin homolog (Drosophila)
GENE ID	23413
ORIGIN SPECIES	Human
ACCESSIONS	BC004856
VECTOR	pDONR223
NCBI SYMBOL	FREQ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaataacacagacaagttggaagcccgaaattgtggagagagctgaccaggaagacactctatccaggaagaaggaggtccagcagttggtaacaaaggtctacataaggactccccagtgggcagctggatgcggcaggcttcagagaatctacaagcaattctccgttcggag
accgccacaagttgtccacattgtttcaacgtcttgaatgaaacaaggacggggcgaaatgagttctccgagttcatccaggcgctgtggctgacctacgggggaacccgtggatgagaagctacggtgggcctcagctctacgacttgacaatgatggcttaccatcaccaggaatgaga
tgctggacattgtggatccattaccagatggtgggggaataccgtggagctccagagaggagagaacactcctgagaagaggtggaccggatcttgcattgatggataagaatgccacgggaagctgacccctgcaggagttccaggaggtcccaaggcagaccggtccattgt
gcagcgtctgccctctcagcagcggcctgtatag

More info: Please contact us for any question: contact@proteogenix.fr