

FAM107B-family with sequence similarity 107, member B Gene

Reference: PTXBC004872

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM107B; C10orf45; family with sequence similarity 107 member B		
For research use only.			

GENE DESCRIPTION	family with sequence similarity 107, member B
GENE ID	83641
ORIGIN SPECIES	Human
ACCESSIONS	BC004872
VECTOR	pDONR223
NCBI SYMBOL	FAM107B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgagccagactacatagaagatgacaatcctgaactcattagggctcagaaactgatcaatcctgtaaaaaccccggaacctcaagattcttcacagagaactcttgaatcaaaaaaggggtcttgcctcagaacaaccagaattgcagaaggtgatgaaaaaagaaa
acgagaccaagtaataaagcagaaggaagaagcacagaagaagaatctgacttggaatagagctattaaacggcagcagaagttggagcagcttgaacttgagaagcagaattgcaagaagagcaagaagaatgcccccagtttgtgaaggtgaaaggcaatctcaggag
aacaggccaagaagtcgccaagcccaggagtcctag