

YIPF4-Yip1 domain family, member 4 Gene

Reference: PTXBC004875

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein YIPF4; FinGER4; Nbla11189; YIP1 family member 4; Yip1 domain family member 4		

For research use only.

GENE DESCRIPTION	Yip1 domain family, member 4
GENE ID	84272
ORIGIN SPECIES	Human
ACCESSIONS	BC004875
VECTOR	pDONR223
NCBI SYMBOL	YIPF4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcctccggcccccggcctatgcccccactaacggggaattcacctttgtctctcagcagacgcggaagatctcagtggttcaatagcatccagatgtcaaattaactctgtggagattttatcaaagaatctacagctactactattctgagacaaagggttatggctg
gcttctggaagttgaagatgatcctgaagatacaagccactcttgaagaattggacattgatctaaaggatatttactacaaaatccgatgtgtttgatgccaatgccatcacttggttttaatagacaagtggtgagagacaatcctgactttggggtcctctggctgtgtttctttctttc
catgatacattataggacagtttaggggtgctcatggattataaccatttggatatgtgttcactaacaattttctactggccagagttctgtgtggaagttgcataatggccaagtcctggagttataggataattcattactctctcattgtaataagccctgtacttttgggtggtggaatcatt
gaagtgggtctacactataaaactgtttgggtgtgtttgggctgcctacagtgctgcttcattgttagtgggtgaagaattcaagacccaaaaagcctctctgattatccaatctttattatacattattttttgtcgttatatactgggtgtgtga