



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SUSD4-sushi domain containing 4 Gene

Reference: PTXBC004888

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO222; sushi domain-containing protein 4; YHGM196; testicular secretory protein Li 55; sushi domain containing 4		

For research use only.

GENE DESCRIPTION	sushi domain containing 4
GENE ID	55061
ORIGIN SPECIES	Human
ACCESSIONS	BC004888
VECTOR	pDONR223
NCBI SYMBOL	SUSD4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]