



Reference: PTXBC004913

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	proline-rich protein PRCC; RCCP1; TPRC; papillary renal cell carcinoma translocation-associated gene protein; papillary renal cell carcinoma (translocation-associated)		

For research use only.

GENE DESCRIPTION	papillary renal cell carcinoma (translocation-associated)
GENE ID	5546
ORIGIN SPECIES	Human
ACCESSIONS	BC004913
VECTOR	pDONR223
NCBI SYMBOL	PRCC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]