

H2AFX-H2A histone family, member X Gene

Reference: PTXBC004915

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2A.X; H2A/X; histone H2AX; H2AX histone; histone H2A.x; H2A histone family member X		

For research use only.

GENE DESCRIPTION	H2A histone family, member X
GENE ID	3014
ORIGIN SPECIES	Human
ACCESSIONS	BC004915
VECTOR	pDONR223
NCBI SYMBOL	H2AFX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgggcccgggaagactggcggcaaggcccgccaaggccaagtcgcgtcgtcgcgcgcgcgcctccagttccagttgggcccgtgtacaccggctgctgcggaaggccactacgccgagcgcgttggcgcggcgccagtgctacgtggcggcagtgctggagtac
ctcaccgctgagatcctggagctggcgggcaatgcggcccgcgacaacaagaagacgcgaatcatccccgccacctgcagctggccatccgcaacgacgaggagctcaacaagctgctggcgcgctgacgatgccaggaggcgtcctgcccaacatccaggccgtgctg
ctgcccaagaagaccagcgccaccgtggggccgaaggcgccctcgggcggcagaaggccaccaggcctccaggagtactaa