

## GPRC5C-G protein-coupled receptor, family C, group 5, member C Gene

Reference: PTXBC004925

|                |                                                                                                                                                                                                                                          |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                      | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                     | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | RAIG-3; RAIG3; G-protein coupled receptor family C group 5 member C; orphan G-protein coupled receptor; retinoic acid responsive gene protein; retinoic acid-induced gene 3 protein; G protein-coupled receptor class C group 5 member C |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | G protein-coupled receptor, family C, group 5, member C                                           |
| <b>GENE ID</b>          | 55890                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC004925                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | GPRC5C                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggccatccacaaagccttgatgtgcctgggactgcctctctctgttccaggggcctgggcccaggccatgtcccacccggctgcagccaaggcctcaacccctgtactacaacctgtgtgaccgctctggggcgtggggcatcgtcctggagccgtggctggggcggg  
cattgtcaccacgtttgtgctcaccatcatcctgtgtggccagcctccctttgtgcaggacacaaagaacggagcctgtctggggacccagctaagaggccgtgtcaccatcacagcgggtacaatgggcagctgctga