

TMEM204-transmembrane protein 204 Gene

Reference: PTXBC004932

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf30; transmembrane protein 204; claudin-like protein 24; claudin-like protein of 24 kDa		

For research use only.

GENE DESCRIPTION	transmembrane protein 204
GENE ID	79652
ORIGIN SPECIES	Human
ACCESSIONS	BC004932
VECTOR	pDONR223
NCBI SYMBOL	TMEM204
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccgtgcagagactcgtggccgcccgtgctggccctgtgtctcactcactcaacaactggcgccctcactcactgggtgtgccagacgtggaggtggcgccagggcagcgtggggctgtggaggtcctgctggctgtggagagaccggggggggcc
gagccctggggccagagccggccaggtggagcgcacatgactgtgagggcgtgggtgtgggtccgagggcagccggtccagaggtccgagggcaccgtcaaacgtcagttcgacatgatgcgcgctgcaacctgtgtggccacggccgctcaccgcagggcagctcacttc
ctcctggggctgtgtggccctgcccctgtgtcaccgacgcccgtgtgggagggagccatggccgtgcattccaactggcgagttttgtcctgggtcgtcgtgactttctacagaattggcccatcaccaacctgtcctgtgtcctgactgaacattggcgctgccttcg
gccacgtggcgccagccatgctcatctggaacattctccacaagggagggactgcatggccccgggtgattgtcagccgctccctgacagcgcgctttgcccgtggcgtggacaatgactacgtggagtcacatgtga