

SCD5-stearoyl-CoA desaturase 5 Gene

Reference: PTXBC004936

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACOD4; FADS4; HSCD5; SCD2; SCD4; stearoyl-CoA desaturase 5; acyl-CoA-desaturase 4; stearoyl-CoA 9-desaturase; stearoyl-CoA desaturase 2; stearoyl-CoA desaturase 4		

For research use only.

GENE DESCRIPTION	stearoyl-CoA desaturase 5
GENE ID	79966
ORIGIN SPECIES	Human
ACCESSIONS	BC004936
VECTOR	pDONR223
NCBI SYMBOL	SCD5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaggcccgccaccgacgccccgaagatcccttctcgacgccaaggaagaaatccgtgccgggctcgaagctctgaggcgccgccccggagaggccaggcgccgccccgagccagacaatcgtctgagggaatgtcgtctgatgagcttgcctcacttg
ggggccgtgtactccctggtgtctatcccaagccaagccactcactctgtctggcctacttctcttctctctgcccgtctgggtgtgacagctggtgccatcgctgtgagccacaggtctaccgggccaagctgcctctgaggatatttctggctgtgccaaactccatggcttt
ccaggcagcagcagaagtggctctgacccagatgcctga