

NOL5A-nucleolar protein 5A (56kDa with KKE/D repeat) Gene

Reference: PTXBC004937

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NOL5A; SCA36; nucleolar protein 56; NOP56 ribonucleoprotein homolog; nucleolar protein 5A (56kDa with KKE/D repeat); NOP56 ribonucleoprotein		

For research use only.

GENE DESCRIPTION	nucleolar protein 5A (56kDa with KKE/D repeat)
GENE ID	10528
ORIGIN SPECIES	Human
ACCESSIONS	BC004937
VECTOR	pDONR223
NCBI SYMBOL	NOL5A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggaagcaatggtcaggcagaggaagcggctgctgagattactaggaagctggagaaacaggagaagaacgcttaagaaggaaagaaacggctggctgcacttgcctcgcgtcttcagaaacagcagtagtactccagaggagtgtgaggagatgagtgaacc
caaaaagaagaaaaagcaaaagccccaggaggttcctcaggagaaatggaatggaagacccatctctcttccaaacccaagaaaaaagaatcttttccaaaggaggagttgatgagtagcgtcttgaagagaccgctggcagcaccagtattccaaaggagaagtctacacc
caaggaggaaacagttatgacctgaggaggcaggccacagaagtggctccaagaaaaagaggaattctccaaaggaggccggtcagcagtgggcctgaaggcggttgcaagagcagctccaagaagaagaaagtccataaagcatccagggaagattag