

## MGC10814-hypothetical protein MGC10814 Gene

Reference: PTXBC004943

|                               |      |                     |                |
|-------------------------------|------|---------------------|----------------|
| Product type:                 | DNA  | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      |      |                     |                |
| <b>For research use only.</b> |      |                     |                |

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | hypothetical protein MGC10814                                                                     |
| <b>GENE ID</b>          | 84757                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC004943                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | MGC10814                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgctgggtcatgggacagcacgtgcaaggacctgtggcaggagccaggagctgtcaggagacagagcaagggaacaggagacttggtgagagacacttccccacgtccagggtga