

PCGF1-polycomb group ring finger 1 Gene

Reference: PTXBC004952

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	2010002K04Rik; NSPC1; RNF3A-2; RNF68; polycomb group RING finger protein 1; RING finger protein 68; nervous system Polycomb-1; polycomb group ring finger 1		

For research use only.

GENE DESCRIPTION	polycomb group ring finger 1
GENE ID	84759
ORIGIN SPECIES	Human
ACCESSIONS	BC004952
VECTOR	pDONR223
NCBI SYMBOL	PCGF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggcttcggaaccagctccagtcagtgataagatggacccgctacggaacgaggagggttcgagtgaaagacaaagactgaatgaacacattgttgcctatgcgcgggctactcgtggatgccaccacatcacagagtgcttcatacttctgcaagagttgtattgtga
agtaacctccaaactagcaagtaactgccccatgtgcaacattaagatccacgagacacagccactgctcaactcaaaactggaccgggtcatgcaggacatcgtgtataagctggtgcctgcttgcagacagtgagagaaacggatcgggaattctaccagtcctcgggtttggaccg
ggtaaccagccctggggaagagccagcactgagcaacctcggcctcccttcagcagcttggaccacttaagccactactatcgtatgatgagcagtgaaacctgtgcctggagcggctgagttctggcaagacaagaataaaagcgtctgcagaacaaagtatgccgatgt
tctgttagagctgaggtacgccatccggagggtcctgtgtcaccgcttgatgctaaacctcagcatgtgcagctccttttgacaatgaagttctccctgatcacatgacaatgaagcagataggtctctccgctgggttcggcaagccatcccttgcctttacaatacagtgatgaaagaga
agaggaggtag