

FAM64A-family with sequence similarity 64, member A Gene

Reference: PTXBC005004

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM64A; CATS; RCS1; CALM interacting protein expressed in thymus and spleen; CALM-interactor expressed in thymus and spleen; regulator of chromosome segregation 1; regulator of chromosome segregation protein 1; family with sequence similarity 64 m		

For research use only.

GENE DESCRIPTION	family with sequence similarity 64, member A
GENE ID	54478
ORIGIN SPECIES	Human
ACCESSIONS	BC005004
VECTOR	pDONR223
NCBI SYMBOL	FAM64A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctctcggggcagacatggggacccctcgtgcgcggagatctctccagcaccaggagcagctgggagacagcaaggagctgcagcctgtggtcagccatcaggagacctctgtagggccctgggggtccctgtgcagacagttcgaaggaggtgccctgagagccgt
caacctcaacctcgcgcagggccctcgtggaacgctggaaacccagagccaggtcagcaggccctccagcgtgcagctcgtcagctaaagagctcttgggtgcccgtgtccagagaatccaggagctcgtccaaagtggcaccagtggtggtggagaccaggtgaagg
ccaggaggcgggaagagaggagcacagaaggcagtgatcccaactcacagcctgagccagaagagcaccggcgtgtctggagccgccccctgcccactcagccgagacccctgggagaaggagacatcaccgctctgtccggatgggtcacatgccaccattacggc
gatcaaggcgggaggtgcttccggagccctactctcaacagagccctctgctctccagcagctctgacagtgacatagagcctgtggggcggaattcagcatctccagaagctgtcccaagagctagatgaagccattatggcggagagagaaacaagccctgtctga
ccgccaaggcttcaactcaaggatgtctatctccctgga