

PHLDA2-pleckstrin homology-like domain, family A, member 2 Gene

Reference: PTXBC005034

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BRW1C; BWR1C; HLDA2; IPL; TSSC3; pleckstrin homology-like domain family A member 2; beckwith-Wiedemann syndrome chromosomal region 1 candidate gene C protein; imprinted in placenta and liver protein; p17-BWR1C; p17-Beckwith-Wiedemann region 1 C; p17-Beckw		

For research use only.

GENE DESCRIPTION	pleckstrin homology-like domain, family A, member 2
GENE ID	7262
ORIGIN SPECIES	Human
ACCESSIONS	BC005034
VECTOR	pDONR223
NCBI SYMBOL	PHLDA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaatccccgacgaggtgctacgcgagggcgagttggagaagcgcagcagcctctccagctatggaagaagaagcgcgggtgctcacctccgaccgcctgagcctgttccccgagccccgcgcgcgccccagagctgcgttccactccatcctcaagtgga
ctgctgtggagcgcacggcgaagtacgtgtacttcacatcgtcaccaccgaccacaaggagatcgacttccgtcgcgcggcgagagctgctggaacgcggccatcgctggtggtctatcgtttccagaaccgcgcgcctgcaggacttgcagccgccaaggaacgacgcg
caccgcccgcaccgcccaggagcctgtggtcgcggcgccgcacccctccgagccctcggagccctccagccatccccgagcccaaacccgcagccatga