

RAB3B-RAB3B, member RAS oncogene family Gene

Reference: PTXBC005035

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB3B, member RAS oncogene family; brain antigen RAB3B; ras-related protein Rab-3B		

For research use only.

GENE DESCRIPTION	RAB3B, member RAS oncogene family
GENE ID	5865
ORIGIN SPECIES	Human
ACCESSIONS	BC005035
VECTOR	pDONR223
NCBI SYMBOL	RAB3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcagtgacagatggtaaaactggagtcgaagatgcctctgaccagaattttgactacatgtttaactgcttatcattggcaacagcagtggtggcaagacctctctctcctcgtatgctgatgacacgttcacccagcctcgttagcaccgtgggcatcgactcaagtgaga
cagttacccgtcacgagaagcgggtgaaactgcagatctgggacacagctgggcaggagcgggtaccggaccatcacacagcctattaccgtggggccatgggcttcattctgatgatgacatcaccaatgaaagatccttcaatgctgtccaagactgggctactagatcaagaccta
ctctgggacaatgcacaagttaattctgggtgggaacaagtgtgacatggaggaagaggggtgttccactgagaaggccagctccttgacagcagcttgggtttgattctttgaagccagtgcaagagagaacatcagtgtaaggcaggcctttgagcgcctgggtgatgccattt
gtgacaagatgtctgattcgtggacacagaccgtcgatctgggtcctccaagaacacgcgtctctcggacacccaccgctgctgcagcagaactgctatgctag