



Reference: PTXBC005044

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NF-E2; p45; transcription factor NF-E2 45 kDa subunit; leucine zipper protein NF-E2; nuclear factor (erythroid-derived 2), 45kD; nuclear factor (erythroid-derived 2), 45kDa; nuclear factor, erythroid-derived 2 45 kDa subunit; p45 NF-E2; nuclear factor, er		

For research use only.

GENE DESCRIPTION	nuclear factor (erythroid-derived 2), 45kDa
GENE ID	4778
ORIGIN SPECIES	Human
ACCESSIONS	BC005044
VECTOR	pDONR223
NCBI SYMBOL	NFE2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]