

DNAJC30-DnaJ (Hsp40) homolog, subfamily C, member 30 Gene

Reference: PTXBC005056

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WBSCR18; dnaJ homolog subfamily C member 30; DnaJ (Hsp40) homolog, subfamily C, member 30; Williams Beuren syndrome chromosome region 18; williams-Beuren syndrome chromosomal region 18 protein; DnaJ heat shock protein family (Hsp40) member C30		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily C, member 30
GENE ID	84277
ORIGIN SPECIES	Human
ACCESSIONS	BC005056
VECTOR	pDONR223
NCBI SYMBOL	DNAJC30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagccatgcgtggcgtggtggcagcggtgttaccttggaggttgcgcagggccgtggccttcacaaaattctgcaccagcctgggcttaagagcagaggacttattccaggcgactgctctattcgcgcacggcgctgtatgatctgctcggcgtccctccacagcca
cgaggcccaaatcaaggcggccttaccgtcagtgcttctctaccacccggaccgcaactccgggagcgcgggagccgccgagcgttcaegcgcacgtctccaggcctacgtggtgctgggcagtgccacccctcgtcgcaagtatgatcgcggcctactcagcgacgagact
ggcgggacgtggcgtccggccctccaggacggcgcacccgacccggcgtcggcgtaccccgccgcccacctctcgaccacagcaggttctcgggcctccccggcgccaaccgcacgatgttcaactttgacgcctctaccaggccactacggggaacaactggagcgg
gaacggcgctgaggggccggggaggcccttcgcaacggcagggatcggtccatgaaaggcctccgctgggagatacccgagacacggctgccattttctcatctttcaatcttcatcatcgcgttttatattaa