

ZNF2-zinc finger protein 2 Gene

Reference: PTXBC005068

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A1-5; ZNF661; Zfp661; zinc finger protein 2; zinc finger protein 2.2; zinc finger protein 661		

For research use only.

GENE DESCRIPTION	zinc finger protein 2
GENE ID	7549
ORIGIN SPECIES	Human
ACCESSIONS	BC005068
VECTOR	pDONR223
NCBI SYMBOL	ZNF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaacagaaagcaagccctcaggggagactcgtgaagaatccctctccgggacaaagccttgcggcgacgggtcagccctgtccagggaaattctcactaaagagagacaccagggaatgcagtgactgtgggaagacctttttgaccactcatccctacccgccatcagag
gactcacactggggagaagccctacgactccgagtggtgggaaagccttcagccagggagcagccctcagcagacatctgatgcacacactggggagagccctacgagtgagtggtgtctcaaaagccttcttgaccgttcgtccctaactgtccatcagcgaattcacactgg
agagaaccccttcagtgcaacgagtggtggaaagcctttttgaccgttcacccctactcgcacaccagagaattcacactggagaagtccttatgaatgtcatcagtggtggaaagcctttagccagaaaagtattctacgcctcagcgaatccacactggcaggaagccttatgagtg
taacgagtgccgggaaagccttctatgtgtctcgtctgaatagacatcagaaagctcatgtggggaccctcgtatcagtgtaacgagtggtggcaagccttcttgaccgtctaccccttacacagcatcagaagatccacactggagacaagcccatgaatgcagcgaatcgggaa
agcctttagccagcgttgcggctcagcggcatcagcgtgtccacacgggagagaagcccttgaatgcactgtgtgtgggaaagtttcagttcaaaatctctgtattcaacatcaacggcgttacgccaacagggaatagactga