

C12orf32-chromosome 12 open reading frame 32 Gene

Reference: PTXBC005106

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf32; HKMT1188; RHINO; RAD9, HUS1, RAD1-interacting nuclear orphan protein 1; Rad9, Rad1, Hus1 interacting nuclear orphan; rAD9, RAD1, HUS1-interacting nuclear orphan protein; RAD9-HUS1-RAD1 interacting nuclear orphan 1		

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 32
GENE ID	83695
ORIGIN SPECIES	Human
ACCESSIONS	BC005106
VECTOR	pDONR223
NCBI SYMBOL	C12orf32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctcccagaaaaaacgcgcagcctcccagaaagcccgctgctgtccaccaaacaccactggaggccccaacacagctgtgcatctacacagcttcccatcactcacactcgacaggtgccagcaagccattgaccacagcaccatcctctgggtatcacctgatt
ttgatacagcagcaggaagctgttcccagcctaccagaaacacaaaacggcgagacactcaagtcgaaaacaccacctccaagttccacatctaacttttgagagtcgcaatctccagttcagagacattgggatcccttaacgagagtgccccagtgaatcagaaaaag
gatgttccagaagacccttagttccagtgctcagtcaccaagctgtgggaacatgtcagtcaggcacttcagagcttaccttatgtgttcattccacctgatatccagacccagagtcacatctgtgaaggagaactcattcccaagatcagaaggaaaacagccttctaagctgcac
tcttcacactggcactcctaatagcccagagcctggagctgttctgttaagacacccccgaggacaagtatggaataaagtcacatggaggagacgacagcactgttgcctacacaggagagagggaagctgagcagaagccaattctgtgaaaagctga